

LUNA™ Fluorescence Calibration Beads

F23102

Storage

4 °C in the dark

Product Description

Appearance	Colorless beads in clear liquid
Bead content	Polystyrene
Excitation/emission	480/520 nm

LUNA™ Fluorescence Calibration Beads are used to calibrate the LUNA™ family of cell counters for fluorescence cell counting.

Directions for Use

Mix LUNA™ Fluorescence Calibration Beads thoroughly before use to ensure a homogenous suspension. Work quickly when using LUNA™ Fluorescence Calibration Beads to ensure the beads do not settle.

Disclaimer

This product is for research use only.
Please consult the material safety data sheet for information regarding hazards and safe handling practices.

Additional information is available on our website at www.logosbio.com.
© Logos Biosystems 2018. All rights reserved.



HEADQUARTERS

FL 2 & 3
28 Simindaero 327beon-gil, Dongan-gu
Anyang-si, Gyeonggi-do 14055
South Korea

Tel: +82 (31) 478-4185

USA

7700 Little River Turnpike STE 207
Annandale, VA 22003
USA

Tel: +1 (703) 622-4660, +1 (703) 942-8867

EUROPE

11B avenue de l'Harmonie
59650 Villeneuve d'Ascq
France

Tel: +33 (0)3 74 09 44 35

www.logosbio.com

VL1801-01