



## Protein Safe Lysing kits

Prevent protein degrading by shutting down unwanted protease activity & ensure your protein stays intact for your experiment



Available for CK-14 & CK28-R lysing kits

### PROTEIN SAFE IS EFFICIENT AGAINST...

- Serine Proteases
- Aminopeptidase B & Leucine Aminopeptidase
- Cysteine Proteases & Trypsin-like Proteases
- Aspartic Proteases

### COMPOSITION

- AEBSF, aprotinin, bestatin, E-64, Leupeptin, Pepstatin A

### PRODUCT REFERENCES


- Protein Safe Hard tissue homogenizing CK28-R:  
**P000972-LYSK0-A.0**
- Protein Safe Soft tissue homogenizing CK14:  
**P000973-LYSK0-A.0**

## PRECELLYS LYSING KITS RANGE



- 5** types of tube volumes  
(0,3mL in well-plate format, 0,5, 2, 7 & 15mL)
- 4** grinding materials  
(ceramic, glass, garnet and stainless steel)
- 30+** lysing kits

## LIFE SCIENCES PRODUCT RANGE



**1500+ SCIENTIFIC DOCUMENTS**

**SIGN UP TO THE APPLICATION CENTER**  
[www.bertin-instruments.com](http://www.bertin-instruments.com)

Bertin Technologies - Head Office  
10 bis Avenue Ampère - Parc d'activités du Pas du Lac  
78180 Montigny-le-Bretonneux FRANCE