



Cell Counters

Fast • Accurate • Affordable

Logos Biosystems

THE STANDARD FOR ACCURATE CELL COUNTING

With a growing emphasis on both inter- and intra-lab reproducibility in research, the seemingly simple act of getting accurate cell counts is now more important than ever.

The automated cell counters from our partner Logos Biosystems are equipped with high-quality optics and sophisticated image analysis software to bring you cell concentration and viability data with incredible speed, accuracy and reliability.

Advanced declustering algorithms distinguish single cells in clusters or clumps. Whether you use disposable precision slides for increased convenience or the budget-friendly LUNA reusable slide, the Logos Biosystems cell counters deliver accurate and reliable counts every single time.

The LUNA-FX7 Automated Cell Counter is the most powerful cell counter to date, with unmatched cell counting accuracy, 1-, 2-, 3- and 8-chamber slide options, state of the art optics, dual fluorescence and brightfield detection, precision autofocus and the option for 21 CFR Part 11/ GMP compliance.

Brightfield or fluorescence

From PBMCs to bacteria

Customizable protocols

Advanced cluster analysis

Convenient workflow

Onboard monitoring

Simple record keeping

www.biocat.com/logos



PBMCs: Peripheral Blood Mononuclear Cells

BioCat GmbH
Technologiepark
Im Neuenheimer Feld 584
D-69120 Heidelberg

Tel.: +49 (0) 6221 71415 16
Fax: +49 (0) 6221 71415 29
E-Mail: info@biocat.com

www.biocat.com

Logos Biosystems Cell Counters



LUNA-III™ Cell Counter

The LUNA-III provides unmatched accuracy with advanced algorithms, precise autofocus, reusable slides, non-toxic dyes, network access, re-analysis options, and a compact, user-friendly design.



LUNA-II™ Cell Counter

Small yet powerful, the LUNA-II is cell counting convenience at its finest. Simply insert your sample and the LUNA-II does the rest: autofocus, adjust light, and count in just 15 seconds.



LUNA-FL™ Dual Fluorescence Cell Counter

Equipped with dual fluorescence and brightfield optics, the LUNA-FL is a powerhouse. The LUNA-FL can assess cell counts, viability, and GFP transfection efficiency without being limited by cell type or size.



LUNA-FX7™ Cell Counter

The LUNA-FX7 offers unmatched cell counting accuracy, a maximum counting volume of 5.1 μ L, all new optics, dual fluorescence and brightfield illumination, precise autofocus and 1/2/3/8-chamber slide options.



LUNA-STEM™ Fluorescence Cell Counter

The LUNA-STEM takes adipose-derived stem cell and stromal vascular fraction samples to count live nucleated cells, dead nucleated cells, and non-nucleated cells with precision and consistency.



QUANTOM Tx™ Microbial Cell Counter

In less than 1 minute, the QUANTOM Tx scans up to 20 fields of view to produce accurate single bacterial cell counts. The advanced software can account for the varying shapes, sizes, and arrangements of different bacteria.

www.biocat.com/logos/compare-counters



BioCat GmbH

Technologiapark
Im Neuenheimer Feld 584
D-69120 Heidelberg

Tel.: +49 (0) 6221 71415 16
Fax: +49 (0) 6221 71415 29
E-Mail: info@biocat.com

www.biocat.com